

# KellMed

**Sterile Packed Disposables for  
Orthopaedic and Trauma Theatres**

---

[www.kellmed.co.uk](http://www.kellmed.co.uk)  
+44 (0) 1305 751821

# Contents

PRODUCT	PAGE	
Wires and Pins	Single Trocar	3
	Double Trocar	5
	Fully Threaded	7
	Threaded Tip	8
	Denham Pins	9
	Trocar / Triangle Round/Round	10
	Cerclage Wire	11
	ACL Guidewire	12
Drills	Quick Fit Connection	13
	Plain Shank	14
	Bone Taps	15
KellMed	Bespoke Sizes	16
	Packing & sterilization	17
	How to order & Delivery	18

# KellMed

**A specialist manufacturer of sterile disposables for orthopaedic and trauma surgery.**

KellMed are a UK family firm, we have been manufacturing and supplying the NHS and private sector since 2002. We provide a range of excellent quality wires, pins, drills and taps at competitive prices. We pride ourselves on impeccable customer service and hold a large range in stock for immediate dispatch. We are also able to offer a bespoke manufacture service for more unusual sizes.

All our products are **supplied sterile, for single use, with a maximum shelf life of five years.** We hold full UKCA and ISO 13485:2016 certification, all products are registered with the MHRA.



**Email**  
office@kellmed.co.uk



**Phone**  
01305 751821



Remember to include your PO number, our product codes and the quantity you require

## Sales Office

Unit 1 Northern Way, Cropmead, Crewkerne, Somerset, TA18 7HJ, UK

Company No: 8040625 | VAT No: GB133087823



# Wires and Pins

## K Wire | Single Trocar

The **single trocar** K wires on this list are the most popular sizes in our range. The majority are available from stock ready for rapid shipment to your stores. If you don't see what you need on this list, let us know. We may well have it and if not, we will be happy to manufacture wires to your specification, and to hold stock to meet your future requirements.

DESCRIPTION	PRODUCT CODE
<b>0.8mm diameter</b>	
0.8mm x 70mm k wire single trocar	319.070s
0.8mm x 80mm k wire single trocar	319.080s
0.8mm x 150mm k wire single trocar	319.150s
<b>0.9mm diameter</b>	
0.9mm x 75mm k wire single trocar	362.075s
0.9mm x 100mm k wire single trocar	362.100s
0.9mm x 125mm k wire single trocar	362.125s
0.9mm x 140mm k wire single trocar	362.140s
0.9mm x 150mm k wire single trocar	362.150s
<b>1.0mm diameter</b>	
1.0mm x 70mm k wire single trocar	320.070s
1.0mm x 100mm k wire single trocar	320.100s
1.0mm x 150mm k wire single trocar	320.150s
<b>1.1mm diameter</b>	
1.1mm x 75mm k wire single trocar	363.075s
1.1mm x 100mm k wire single trocar	363.100s
1.1mm x 125mm k wire single trocar	363.125s
1.1mm x 150mm k wire single trocar	363.150s
1.1mm x 200mm k wire single trocar	363.200s
<b>1.2mm diameter</b>	
1.2mm x 150mm k wire single trocar	321.150s
1.2mm x 250mm k wire single trocar	321.250s
1.2mm x 300mm k wire single trocar	321.300s
1.2mm x 305mm k wire single trocar	321.305s
<b>1.25mm diameter</b>	
1.25mm x 125mm k wire single trocar	346.125s
1.25mm x 150mm k wire single trocar	346.150s
1.25mm x 280mm k wire single trocar	346.280s
<b>1.4mm diameter</b>	
1.4mm x 100mm k wire single trocar	322.100s
1.4mm x 125mm k wire single trocar	322.125s
1.4mm x 150mm k wire single trocar	322.150s

DESCRIPTION	PRODUCT CODE
<b>1.5mm diameter</b>	
1.5mm x 125mm k wire single trocar	364.125s
1.5mm x 150mm k wire single trocar	364.150s
1.5mm x 200mm k wire single trocar	364.200s
1.5mm x 230mm k wire single trocar	364.230s
1.5mm x 250mm k wire single trocar	364.250s
1.5mm x 300mm k wire single trocar	364.300s
1.5mm x 310mm k wire single trocar	364.310s
<b>1.6mm diameter</b>	
1.6mm x 75mm k wire single trocar	323.075s
1.6mm x 100mm k wire single trocar	323.100s
1.6mm x 125mm k wire single trocar	323.125s
1.6mm x 150mm k wire single trocar	323.150s
1.6mm x 200mm k wire single trocar	323.200s
1.6mm x 230mm k wire single trocar	323.230s
1.6mm x 280mm k wire single trocar	323.280s
<b>1.8mm diameter</b>	
1.8mm x 150mm k wire single trocar	325.150s
<b>2.0mm diameter</b>	
2.0mm x 75mm k wire single trocar	324.075s
2.0mm x 100mm k wire single trocar	324.100s
2.0mm x 125mm k wire single trocar	324.125s
2.0mm x 150mm k wire single trocar	324.150s
2.0mm x 200mm k wire single trocar	324.200s
2.0mm x 230mm k wire single trocar	324.230s
2.0mm x 280mm k wire single trocar	324.280s
2.0mm x 310mm k wire single trocar	324.310s
<b>2.4mm diameter</b>	
2.4mm x 200mm steinmann pin single trocar	334.200s
<b>2.5mm diameter</b>	
2.5mm x 150mm k wire single trocar	335.150s
2.5mm x 175mm k wire single trocar	335.175s
2.5mm x 200mm k wire single trocar	335.200s
2.5mm x 250mm k wire single trocar	335.250s
<b>2.8mm diameter</b>	
2.8mm x 230mm steinmann pin single trocar	309.230s
2.8mm x 250mm steinmann pin single trocar	309.250s
<b>3.2mm diameter</b>	
3.2mm x 200mm k wire single trocar	312.200s
3.2mm x 230mm steinmann pin single trocar	312.230s
3.2mm x 250mm steinmann pin single trocar	312.250s

### All our wires are:

- > Wires are individually double pouched in tough Tyvek® and **sold in packs of 10.**
- > Manufactured from 316L implant grade stainless steel.
- > Each outer pouch carries three self-adhesive traceability labels featuring product code, LOT number, and expiry date.
- > All our wires and pins are UKCA marked and MHRA registered.

# Wires and Pins

## K Wire | Double Trocar

The **double trocar** K wires on this list are the most popular sizes in our range. The majority are available from stock ready for rapid shipment to your stores. If you don't see what you need on this list, let us know. We may well have it and if not, we will be happy to manufacture wires to your specification, and to hold stock to meet your future requirements.

DESCRIPTION	PRODUCT CODE
<b>0.6mm diameter</b>	
0.6mm x 100mm k wire double trocar	232.100s
0.6mm x 150mm k wire double trocar	232.150s
<b>0.8mm diameter</b>	
0.8mm x 80mm k wire double trocar	338.080s
0.8mm x 100mm k wire double trocar	338.100s
0.8mm x 150mm k wire double trocar	338.150s
<b>0.9mm diameter</b>	
0.9mm x 75mm k wire double trocar	370.075s
0.9mm x 100mm k wire double trocar	370.100s
0.9mm x 125mm k wire double trocar	370.125s
0.9mm x 130mm k wire double trocar	370.130s
0.9mm x 150mm k wire double trocar	370.150s
<b>1.0mm diameter</b>	
1.0mm x 70mm k wire double trocar	339.070s
1.0mm x 100mm k wire double trocar	339.100s
1.0mm x 125mm k wire double trocar	339.125s
1.0mm x 150mm k wire double trocar	339.150s
<b>1.1mm diameter</b>	
1.1mm x 75mm k wire double trocar	371.075s
1.1mm x 100mm k wire double trocar	371.100s
1.1mm x 125mm k wire double trocar	371.125s
1.1mm x 150mm k wire double trocar	371.150s
1.1mm x 200mm k wire double trocar	371.200s
1.1mm x 300mm k wire double trocar	371.300s
<b>1.2mm diameter</b>	
1.2mm x 100mm k wire double trocar	340.100s
1.2mm x 125mm k wire double trocar	340.125s
1.2mm x 150mm k wire double trocar	340.150s
1.2mm x 200mm k wire double trocar	340.200s
1.2mm x 230mm k wire double trocar	340.230s

DESCRIPTION	PRODUCT CODE
<b>1.25mm diameter</b>	
1.25mm x 75mm k wire double trocar	345.075s
1.25mm x 100mm k wire double trocar	345.100s
1.25mm x 125mm k wire double trocar	345.125s
1.25mm x 150mm k wire double trocar	345.150s
<b>1.4mm diameter</b>	
1.4mm x 75mm k wire double trocar	341.075s
1.4mm x 100mm k wire double Trocar	341.100s
1.4mm x 125mm k wire double trocar	341.125s
1.4mm x 150mm k wire double trocar	341.150s
<b>1.5mm diameter</b>	
1.5mm x 100mm k wire double trocar	372.100s
1.5mm x 150mm k wire double trocar	372.150s
<b>1.6mm diameter</b>	
1.6mm x 75mm k wire double trocar	342.075s
1.6mm x 100mm k wire double trocar	342.100s
1.6mm x 125mm k wire double trocar	342.125s
1.6mm x 140mm k wire double trocar	342.140s
1.6mm x 150mm k wire double trocar	342.150s
1.6mm x 175mm k wire double trocar	342.175s
1.6mm x 200mm k wire double trocar	342.200s
1.6mm x 250mm k wire double trocar	342.250s
1.6mm x 280mm k wire double trocar	342.280s
<b>1.8mm diameter</b>	
1.8mm x 125mm k wire double trocar	343.125s
1.8mm x 150mm k wire double trocar	343.150s
<b>2.0mm diameter</b>	
2.0mm x 75mm k wire double trocar	344.075s
2.0mm x 100mm k wire double trocar	344.100s
2.0mm x 125mm k wire double trocar	344.125s
2.0mm x 150mm k wire double trocar	344.150s
2.0mm x 200mm k wire double trocar	344.200s
2.0mm x 230mm k wire double trocar	344.230s
<b>2.5mm diameter</b>	
2.5mm x 150mm k wire double trocar	337.150s

**All our wires are:**

- > Wires are individually double pouched in tough Tyvek® and **sold in packs of 10.**
- > Manufactured from 316L implant grade stainless steel.
- > Each outer pouch carries three self-adhesive traceability labels featuring product code, LOT number, and expiry date.
- > All our wires and pins are UKCA marked and MHRA registered.

# Wires and Pins

## Threaded Wires

We carry a range of threaded wires, partially with a 15mm tip, fully threaded – single and double trocar, and Denham Pins with a centrally threaded section.

The guide wires on this list are the most popular sizes in our range. They are all held in stock ready for rapid shipment to your stores. If you don't see what you need on this list, let us know, we may have it in stock. If not, we will be happy to manufacture wires to your specification, and to hold stock to meet your future requirements.

## Wires and Pins Fully Threaded

DESCRIPTION	PRODUCT CODE
<b>1.5mm diameter</b>	
1.5mm x 100mm k wire double trocar fully threaded	316.100s
<b>1.6mm diameter</b>	
1.6mm x 100mm k wire double trocar fully threaded	311.100s
1.6mm x 125mm k wire double trocar fully threaded	311.125s
<b>2.0mm diameter</b>	
2.0mm x 230mm steinmann pin single trocar fully threaded	381.230s
<b>2.4mm diameter</b>	
2.4mm x 230mm steinmann pin single trocar fully threaded	317.230s
<b>2.5mm diameter</b>	
2.5mm x 230mm steinmann pin single trocar fully threaded	382.230s

DESCRIPTION	PRODUCT CODE
<b>2.8mm diameter</b>	
2.8mm x 230mm steinmann pin single trocar fully threaded	318.230s
<b>3.0mm diameter</b>	
3.0mm x 200mm steinmann pin single trocar fully threaded	383.200s
3.0mm x 230mm steinmann pin single trocar fully threaded	383.230s
<b>3.2mm diameter</b>	
3.2mm x 230mm steinmann pin single trocar fully threaded	384.230s
<b>3.5mm diameter</b>	
3.5mm x 230mm steinmann pin single trocar fully threaded	385.230s
<b>4.0mm diameter</b>	
4.0mm x 200mm steinmann pin single trocar fully threaded	386.200s
4.0mm x 230mm steinmann pin single trocar fully threaded	386.230s

### All our wires are:

- > Wires are individually double pouched in tough Tyvek® and **sold in packs of 10.**
- > Manufactured from 316L implant grade stainless steel.
- > Each outer pouch carries three self-adhesive traceability labels featuring product code, LOT number, and expiry date.
- > All our wires and pins are UKCA marked and MHRA registered.

## Threaded Tip

DESCRIPTION	PRODUCT CODE
<b>1.25mm diameter</b>	
1.25mm x 150mm k wire single trocar threaded tip	347.150s
<b>1.3mm diameter</b>	
1.3mm x 150mm guide wire single trocar with threaded tip	373.150s
<b>1.4mm diameter</b>	
1.4mm x 150mm guide wire single trocar with threaded tip	375.150s
<b>1.6mm diameter</b>	
1.6mm x 150mm guide wire single trocar threaded tip	348.150s
<b>2.0mm diameter</b>	
2.0mm x 150mm k wire single trocar threaded tip	359.150s
2.0mm x 230mm guide wire single trocar with threaded tip	359.230s
2.0mm x 280mm guide wire single trocar with threaded tip	359.280s
<b>2.4mm diameter</b>	
2.4mm x 229mm k wire threaded trocar tip	356.229s
<b>2.5mm diameter</b>	
2.5mm x 230mm k wire single trocar threaded tip	353.230s
2.5mm x 260mm k wire single trocar threaded tip	353.260s
<b>3.0mm diameter</b>	
3.0mm x 150mm k wire single trocar with threaded tip	351.150s
<b>3.2mm diameter</b>	
3.2mm x 300mm guide wire single trocar threaded tip	376.300s
3.2mm x 400mm guide wire single trocar threaded tip	376.400s

### All our wires are:

- > Wires are individually double pouched in tough Tyvek® and **sold in packs of 10.**
- > Manufactured from 316L implant grade stainless steel.
- > Each outer pouch carries three self-adhesive traceability labels featuring product code, LOT number, and expiry date.
- > All our wires and pins are UKCA marked and MHRA registered.



# Wires and Pins

## Denham Pins

DESCRIPTION	PRODUCT CODE
<b>4.0mm diameter</b>	
4.0mm x 200mm denham pin	390.200s
4.0mm x 230mm denham pin	390.230s
<b>4.8mm diameter</b>	
4.8mm x 230mm denham pin	391.230s
<b>5.0mm diameter</b>	
5.0mm x 200mm denham pin	397.200s
5.0mm x 230mm denham pin	397.230s

### All our wires are:

- > Wires are individually double pouched in tough Tyvek® and **sold in packs of 10.**
- > Each outer pouch carries three self-adhesive traceability labels featuring product code, LOT number, and expiry date.
- > All our wires and pins are UKCA marked and MHRA registered.

# Wires and Pins

## Steinmann Pin Trocar/Triangle

The Steinmann Pins on this list are the most popular sizes in our range. If you don't see what you need on this list, let us know, we may have it in stock. If not, we will be happy to manufacture to your specification.

DESCRIPTION	PRODUCT CODE
<b>2.8mm diameter</b>	
2.8mm x 200mm steinmann pin trocar/triangle	308.200s
<b>3.5mm diameter</b>	
3.5mm x 200mm steinmann pin trocar/triangle	313.200s
<b>4.0mm diameter</b>	
4.0mm x 225mm steinmann pin trocar/ triangle	306.225s
4.0mm x 254mm steinmann pin trocar/ triangle	306.254s
4.0mm x 150mm steinmann pin trocar/ triangle	306.150s
4.0mm x 125mm steinmann pin trocar/ triangle	306.125s
4.0mm x 175mm steinmann pin trocar/ triangle	306.175s
4.0mm x 200mm steinmann pin trocar/ triangle	306.200s

# Wires and Pins

## Round/Round

DESCRIPTION	PRODUCT CODE
<b>1.6mm diameter</b>	
1.6mm x 150mm k wire round /round	404.150s
<b>2.0mm diameter</b>	
2.0mm x 150mm k wire round/round	401.150s
2.0mm x 200mm k wire round/round	401.200s

### All our wires are:

- > Wires are individually double pouched in tough Tyvek® and **sold in packs of 10.**
- > Each outer pouch carries three self-adhesive traceability labels featuring product code, LOT number, and expiry date.
- > All our wires and pins are UKCA marked and MHRA registered.

# Wires and Pins

## Cerclage Wire

All our **Cerclage Wire** is sold in ready cut lengths, individually pouched, sterile packed with full traceability.

DESCRIPTION	PRODUCT CODE
<b>0.4mm diameter</b>	
0.4mm x 400mm cerclage bone wire	720.400s
<b>0.6mm diameter</b>	
0.6mm x 400mm cerclage bone wire	721.400s
0.6mm x 500mm cerclage bone wire	721.500s
<b>0.7mm diameter</b>	
0.7mm x 400mm cerclage bone wire	761.400s
0.7mm x 500mm cerclage bone wire	761.500s
<b>0.8mm diameter</b>	
0.8mm x 400mm cerclage bone wire	730.400s
0.8mm x 500mm cerclage bone wire	730.500s
<b>0.9mm diameter</b>	
0.9mm x 400mm cerclage bone wire	731.400s
<b>1.0mm diameter</b>	
1.0mm x 500mm cerclage bone wire	735.500s
1.0mm x 600mm cerclage bone wire	735.600s

DESCRIPTION	PRODUCT CODE
<b>1.2mm diameter</b>	
1.2mm x 300mm cerclage bone wire	750.300s
1.2mm x 400mm cerclage bone wire	750.400s
1.2mm x 500mm cerclage bone wire	750.500s
1.2mm x 600mm cerclage bone wire	750.600s
looped cerclage wire 1.2 x 900mm	750.090s
<b>1.5mm diameter</b>	
1.5mm x 400mm cerclage bone wire	741.400s
1.5mm x 500mm cerclage bone wire	741.500s
<b>2.0mm diameter</b>	
2.0mm x 300mm cerclage bone wire	790.300s
2.0mm x 400mm cerclage bone wire	790.400s

# Wires and Pins

## ACL Guide Wires, Suture Eye with Drill Tip

DESCRIPTION	PRODUCT CODE
ACL Guide Wire with Suture Eye Drill Tip 2.4mm x 380mm	610.380s

### All our wires are:

- > Individually double pouched in tough Tyvek® and **sold in packs of 10.**
- > Manufactured from 316L implant grade stainless steel
- > Each outer pouch carries three self-adhesive traceability labels featuring product code, LOT number, and expiry date.
- > UKCA marked and MHRA registered.

### All our wires are:

- > Wires are individually double pouched in tough Tyvek® and **sold in packs of 10.**
- > Each outer pouch carries three self-adhesive traceability labels featuring product code, LOT number, and expiry date.
- > All our wires and pins are UKCA marked and MHRA registered.



**Email**  
office@kellmed.co.uk



**Phone**  
01305 751821



Remember to include your PO number, our product codes and the quantity you require

# Drills

## Quick Fit Connection



DESCRIPTION	PRODUCT CODE
<b>1.5mm diameter</b>	
1.5mm x 85mm drill bit quick fit	913.085s
<b>1.8mm diameter</b>	
1.8mm x 100mm drill bit quick fit	920.100s
<b>2.0mm diameter</b>	
2.0mm x 100mm drill bit quick fit	914.100s
2.0mm x 125mm drill bit quick fit	914.125s
2.0mm x 75mm drill bit quick fit	914.075s
<b>2.4mm diameter</b>	
2.4mm x 100mm drill bit quick fit	921.100s
2.5mm x 110mm drill bit quick fit	915.110s
<b>2.5mm diameter</b>	
2.5mm x 125mm drill bit quick fit	915.125s
2.5mm x 170mm drill bit quick fit	915.170s
<b>2.7mm diameter</b>	
2.7mm x 100mm drill bit quick fit	916.100s
2.7mm x 125mm drill bit quick fit	916.125s
<b>3.0mm diameter</b>	
3.0mm x 190mm drill bit quick fit	983.190s
<b>3.2mm diameter</b>	
3.2mm x 125mm drill bit quick fit	917.125s
3.2mm x 145mm drill bit quick fit	917.145s
3.2mm x 195mm drill bit quick fit	917.195s

DESCRIPTION	PRODUCT CODE
<b>3.5mm diameter</b>	
3.5mm x 110mm drill bit quick fit	918.110s
3.5mm x 145mm drill bit quick fit	918.145s
3.5mm x 195mm drill bit quick fit	918.195s
<b>4.5mm diameter</b>	
4.5mm x 145mm drill bit quick fit	919.145s
4.5mm x 195mm drill bit quick fit	919.195s

# Drills

## Plain Shank



DESCRIPTION	PRODUCT CODE
<b>1.5mm diameter</b>	
1.5mm x 70mm drill bit plain shank	931.070s
<b>2.0mm diameter</b>	
2.0mm x 85mm drill bit plain shank	932.085s
<b>2.4mm diameter</b>	
2.4mm x 125mm drill bit plain shank	980.125s
<b>2.5mm diameter</b>	
2.5mm x 95mm drill bit plain shank	933.095s
2.5mm x 180mm drill bit plain shank	933.180s
<b>2.7mm diameter</b>	
2.7mm x 100mm drill bit plain shank	934.100s
<b>2.8mm diameter</b>	
2.8mm x 125mm drill bit plain shank	981.125s
<b>4.0mm diameter</b>	
4.0mm x 120mm drill bit plain shank	939.120s
4.0mm x 180mm drill bit plain shank	939.180s
<b>4.5mm diameter</b>	
4.5mm x 180mm drill bit plain shank	937.180s
<b>6.5mm diameter</b>	
6.5mm x 150mm drill bit plain shank	938.150s

### All our drills are:

- > Drills are individually double pouched in tough Tyvek® and **sold in packs of 5.**
- > Each outer pouch carries three self-adhesive traceability labels featuring product code, LOT number, and expiry date.
- > All our drills are UKCA marked and MHRA registered.

### All our drills are:

- > Drills are individually double pouched in tough Tyvek® and **sold in packs of 5.**
- > Each outer pouch carries three self-adhesive traceability labels featuring product code, LOT number, and expiry date.
- > All our drills UKCA marked and MHRA registered.

# Drills

## Bone Taps

DESCRIPTION	PRODUCT CODE
<b>3.5mm diameter</b>	
cortical tap 3.5mm x 110mm quick fit	941.110s
<b>4.0mm diameter</b>	
cancellous tap 4.0mm x 110mm quick fit	945.110s
<b>4.5mm diameter</b>	
cortical tap 4.5mm x 125mm quick fit	942.125s



### All our drills are:

- > Taps are individually double pouched in tough Tyvek® and **sold in packs of 5.**
- > Each outer pouch carries three self-adhesive traceability labels featuring product code, LOT number, and expiry date.
- > All our taps are UKCA marked and MHRA registered.



Email

office@kellmed.co.uk



Phone

01305 751821



Remember to include your PO number, our product codes and the quantity you require

# Bespoke Sizes


- We are able to offer a custom manufacture service for sizes that are not part of our standard ranges. If you can't see what you need in the catalogue please get in touch to discuss your requirements.
- We have a **low minimum order quantity** and goods will be manufactured, packed and sterilized within 14-16 weeks.
- If required, they can be held in stock thereafter for future orders.
- Please telephone or email us with your requirements and we will provide a quotation.



office@kellmed.co.uk




01305 751821




# How to order




  
**Email**  
 office@kellmed.co.uk



  
**Phone**  
 01305 751821



  
 Remember to include your PO number, our product codes and the quantity you require

## Sales Office

Unit 1 Northern Way, Cropmead,  
 Crewkerne, Somerset, TA18 7HJ, UK

Company No: 8040625  
 VAT No: GB133087823



# Delivery

**Most sizes in this catalogue are held in stock for immediate dispatch.**

Orders received before midday can be dispatched the same day on a 24hr courier service if required.

Orders below £250.00 incur a minimal delivery charge, order with a value above £250.00 are delivered free of charge.

All prices quoted ex VAT



# Packing and sterilization

All our products are packed in inner and outer pouches of 75gsm tough Tyvek with a clear polyethylene front, heat sealed and fully compliant to recognised standards.

All products are sterilized with a 25 kGy gamma radiation dose and have a five-year shelf life post sterilization. Each individual pouch carries three self-adhesive traceability labels featuring product code, LOT number and expiry date.



# KellMed

**Sterile Packed Disposables for  
Orthopaedic and Trauma Theatres**

[www.kellmed.co.uk](http://www.kellmed.co.uk)

+44 (0) 1305 751821

**Sales Office**

Unit 1 Northern Way, Cropmead,  
Crewkerne, Somerset, TA18 7HJ

Company No: 8040625

VAT No: GB133087823

